

IMPORTANT FOR FUTURE REFERENCE	
Please complete this information and retain this manual for the life of the equipment:	
Model #:	_____
Serial #:	_____
Date Purchased:	_____

Installation & Operation Manual

Sterilmatic - Digital Electric Sterilizer STM-ED & STM-EDX



STM-ED

⚠ WARNING

Improper installation, adjustment, alteration, service or maintenance can cause property damage, injury or death. Read the installation, operating and maintenance instructions thoroughly before installing or servicing this equipment.

CROWN FOOD SERVICE EQUIPMENT

 **MIDDLEBY** A Middleby Company

70 Oakdale Road, Downsview (Toronto) Ontario, Canada, M3N 1V9

Telephone: 919-762-1000

www.crownsteamgroup.com

Your Service Agency's Address:

Model

Serial number

Kettle installed by

Installation checked by

SAFETY PRECAUTIONS

Before installing and operating this equipment, be sure everyone involved in its operation is fully trained and aware of precautions. Accidents and problems can be caused by failure to follow fundamental rules and precautions.

The following symbols, found throughout this manual, alert you to potentially dangerous conditions to the operator, service personnel, or to the equipment.



This symbol warns of immediate hazards that will result in severe injury or death.



This symbol refers to a potential hazard or unsafe practice that could result in injury or death.



This symbol refers to a potential hazard or unsafe practice that could result in injury, product damage, or property damage.



This symbol refers to information that needs special attention or must be fully understood, even though not dangerous.

IMPORTANT NOTES FOR INSTALLATION AND OPERATION

WARNING

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury or death.

WARNING

FOR YOUR SAFETY:

Do not store or use gasoline or other flammable vapors or liquids in the vicinity of this or any other appliance.

WARNING

Improper installation, operation, adjustment, alteration, service or maintenance can cause property damage, injury or death. Read the installation, operating and maintenance instructions thoroughly before installing, operating or servicing this equipment.

NOTICE

This product is intended for commercial use only. NOT FOR HOUSEHOLD USE.

NOTICE

This manual should be retained for future reference.

IMPORTANT

The information contained in this manual is important for the proper installation, use, and maintenance of this kettle. Adherence to these procedures and instructions will result in satisfactory baking results and long, trouble free service. Please read this manual carefully and retain it for future reference.


ERRORS: Descriptive, typographic or pictorial errors are subject to correction. Specifications are subject to change without notice.

RETAIN THIS MANUAL FOR FUTURE REFERENCE.

Table of Contents

Important Notes For Installation and Operation	3
Service Connections	5
Introduction	6
Installation	8
Operation	13
Maintenance	21
Warranty	22

SERVICE CONNECTIONS



SERVICE CONNECTIONS

SYMBOL	DESCRIPTION	NOTES
A	Drain - 1/2" (13mm) N.P.T. of 3/8" (10mm) OD copper (see note 1)	1. An air break must be provided if a unit drain line is run. 2. Vent exhaust to atmosphere. B1 is actual connection, but must exit casing at B. 3. Heat of Rejection (Heat Loss Into the Room) 6,000 BTU'S. 4. For use at altitudes up to 6500ft (200m).
B	Steam Exhaust Connection - 5/8" (16mm) (see note 2)	
C	Electrical Connection - (*see electrical specifications)	
D	Power Supply	

The manufacturer reserves the right to modify materials and specifications without notice.

ELECTRICAL CHARACTERISTICS

Model	kW	Hz	1 Phase			3 Phase			
			208V	220V	240V	208V	240V	220/380V	240/415V
STM-ED	9.3	60	45A	-	-	26A	-	-	-
	12.4	60	-	-	52A	-	30A	-	-
STM-EDX	10.4	50	-	48A	-	-	-	16A	-
	12.4	50	-	-	52A	-	-	-	18A

Unit must be grounded. Main supply voltage fluctuations are not to exceed $\pm 10\%$ nominal supply voltage. Clearance from combustible and non-combustible construction should be a minimum of 4" in the back and 2" on the sides.

WATER SUPPLY AND DRAIN SPECIFICATIONS

Good quality water feed is the responsibility of the owner. Water quality must be within the following general guidelines.

TDS: 40-125 ppm Hardness: 35-100 ppm pH: 7.0 - 8.5 Silica: <13 ppm Chlorides: <25 ppm Chlorine: <0.2 ppm Chloramine: <0.2 ppm

The best defense against poor water quality is a water treatment system designed to meet your water quality conditions.

Appliance to be installed with backflow protection according to federal, state or local codes.

IMPORTANT: Exhaust line must be vented to the outside to eliminate the exhaust steam and the accompanying noise from entering the room. Use 1/2" (13mm) copper tubing or suitable alternate. The overall height and length of the line should not rise more than 4' (1.2 meters) above the unit and exceed 15' (4.5 meters) with a minimum of bends. The line should slope downward after leaving the sterilizer in order to ensure condensate drainage.

IMPORTANT: Failure to comply with this outline will affect the sterilization process.

When an exhaust condenser is supplied; the following services must be provided: 1/2" (13mm) cold water: 1" (25mm) waste: 115V electrical line.

DISCLAIMER

Terry System Cartridge Changes / Installation – 2-3 gallons of water MUST be purged at each cartridge change or new installation prior to water supply being fed to the steamer. Failure to do so can result in component damage within the steamer which is not covered under warranty. For additional guidance on proper installation, refer to install documentation provided with each Terry System and Replacement Cartridge Set."

As continued product improvement is a policy of Market Forge, specifications are subject to change without notice.

INTRODUCTION

PRODUCT DESCRIPTION

The Market Forge Sterilmatic Sterilizer with Digital Controller (model STM-ED) is a compact, automatic, low cost steam pressure sterilizer (autoclave).

The sterilizing cylinder is a 3/16" (4.8mm) thick wall, welded aluminum. The exterior is made of polished stainless steel. Interior dimensions of 16" (406mm) diameter and 26" (660mm) long with a cubic content of 5,220 cubic

inches (0.085 cubic meters) and has a door opening of 13-1/2" (343mm) wide and 11" (279mm) high. The sterilizing compartment has a pan capacity of:

- (3) 12" x 20" x 2 1/2" (305mm x 508mm x 64mm) or,
- (2) 12" x 20" x 4" (305mm x 508mm x 102mm) or,
- (1) 12" x 20" x 6" (305mm x 508mm x 152mm)

The sterilizer door is a self-sealing type that cannot be opened until the steam pressure is completely exhausted from the chamber. The door is 12 gauge stainless steel and removable for cleaning without tools. The door gasket is one-piece molded, also replaceable without tools or cement.

The sterilizing cycle is fully automatic with the time, temperature and venting controlled by the microprocessor based, digital controller.

The sterilizing temperature can be set anywhere in a range from 225°F (107°C) and 250°F (121°C). There is an on-board data logger/printer. The data logger records the time, temperature and pressure for each sterilization cycle. This data can be stored for future printing or printed out following each cycle.

OPERATING ENVIRONMENTAL CONDITIONS

This unit is designed for commercial use and to be safe at least under the following conditions:

- For indoor use only.
- For use at altitudes up to 6500ft (2000m)
- For use at temperatures from 41°F (5°C) to 104°F (40°C).
- Maximum relative humidity 80% for temperatures up to 88°F (31°C) decreasing linearly to 50% relative humidity at 104°F (40°C).
- Main supply voltage fluctuations not to exceed $\pm 10\%$ of nominal voltage.
- Transient overvoltages according to Installation Categories II (in accordance with IEC 664).
- Pollution Degree 2 (in accordance with IEC 664).

SERVICE & TECHNICAL INFORMATION CONTACT

WARNING

This unit should be serviced by qualified service personnel only.

Your Sterilmatic Sterilizer has been developed to answer the need for a compact, automatic, low-cost steam pressure sterilizer. The following instructions cover installation.

Should service be required, it is readily available by contacting our authorized service agency located nearest to you. The name of your local service company can be obtained on our website, www.mfii.com.

INTENDED USE - STERILIZATION CYCLE

This unit is intended to be operated intermittently. After a pre-heat cycle, the longest period of sterilization (heating) should be a maximum of 60 minutes. The digital timer allows up to 99 minutes, but it should be kept at no more than 60 minutes. After each use the unit should be opened for removal and reloading of product. The water level should be checked after each use and refilled when necessary.

WATER CONDITIONS

When sterilizing culture mediums that generate sulfide gas or chlorine gas, the inside of the chamber must be cleaned and rinsed thoroughly without fail.

Market Forge from time to time is asked the question about using distilled or deionized water for use with our Sterilizer models STM-ED and STM-EDX. We are always asked why these water choices are not allowed for use with our units and what would be recommended. To address this situation, we have compiled the following as a means of satisfying these questions:

1. We have found that the use of distilled or deionized water will aggressively attack the pure coat of Aluminum Alclad, which protects the bottom surface from oxidizing and then eventually pitting (reference: Operating and Maintenance Instructions).
2. In addition pitting can also be caused by several other external environmental factors. Few examples are as follows. These conditions have been highlighted in our documentation.
 - Grains of hardness in the water supply should be as follows (.25 to 2).
 - A pH imbalance in the water supply can greatly affect the life to the aluminum cylinder. The pH range that would be recommended is between 7.0-8.5.
 - The lack of a positive electrical ground can cause an electrolytic reaction that will accelerate pitting.
 - Another contribution to accelerate pitting is the type of cleaning solutions used or the abrasive scrubbing pads. If a low pH is present with the detergents being used or an abrasive pad, the protective Alclad coating will be removed during the cleaning process.
 - Spillage of media being sterilized can also contribute to the accelerated pitting if it is corrosive.
 - CHLORINE LEVEL \leq 1 PPM.

IMPORTANT

Market Forge will not be responsible for damage resulting from the use of hard or corrosive water, from failure to drain the unit daily, or from inadequate cleaning procedures.

INSTALLATION INSTRUCTIONS

INSTALLATION

Before you begin the installation, check the sterilizer for any transport damages. If there are signs of damage, inform your supplier and freight company immediately.

We recommend using a pallet trolley for lifting the sterilizer. Each model weights 150 lbs. Use trolley to move sterilizer to the place of desired installation, raise the sterilizer slightly higher than counter or stand. Set sterilizer on counter, using the 6" (152mm) legs provided or assembles the optional stainless steel stand with undershelf (follow provided L-1132 instructions). Clearance from combustible and non-combustible construction should be a minimum of 4" in the back and 2" on the sides. If your Sterilmatic includes a water-cooled exhaust condenser, we recommend the use of the Sterilmatic stand, part number 95-6060. First, level unit in place, then adjust rear legs to pitch the unit forward 1/4" (6mm) to insure positive drainage of the cylinder.

ELECTRICAL

Connect to proper electrical supply box and disconnect switch as shown on one of the following schematic diagrams - 208 or 240 volts, single or three phase. Connection is made from the rear of the unit, through the conduit to the terminal box located at the front of the unit.

Electric supply connection for STM-ED

Connect to proper electrical supply as indicated on nameplate on top of unit. The power supply cord is brought in from the rear of the unit, through the conduit and the connection is made at the terminal box located at the front of the unit.

Electric supply connection for STM-EDX (export models)

Connect to proper electrical supply as indicated on nameplate on top of unit. Connection is made from the rear of the unit, through the conduit to the terminal box located at the front of the unit. All control circuits are 220 volts.

In order to accomplish this, a current-carrying grounded neutral must be provided.

Thus, a three phase system must be 4-wires. Most electrical codes require, and we recommend, that a separate switch be located within sight of the sterilizer.

OUTSIDE VENTING

Connect 1/2" (13mm) nominal tubing exhaust to outside vent connection located on top of unit, within the control housing.

IMPORTANT

Exhaust line must be vented to the outside to eliminate the exhausted steam and the accompanying noise from entering the room.

Use 1/2" (13mm) copper tubing or suitable alternate. The overall height and length of the line should not rise more than 4 feet (1.2 meters) above the unit and exceed 15 feet (4.5 meters) with a minimum of bends. The line should slope downward after leaving the sterilizer in order to insure condensate drainage.

WATER-COOLED EXHAUST CONDENSER

If outside venting is not possible, an optional water-cooled condenser is available for connection to an open drain. If required order part no. 95-0436 kit.

TRAY SUPPORTS

Install side tray supports. Tray supports are attached by means of key-hole clearance slots which are slipped over studs located on the sides of the Sterilmatic chamber.

BAFFLE INSTALLATION

To insure maximum drying of packs, a baffle is supplied with your Sterilmatic. Place perforated splash baffle in bottom of the sterilizing chamber. Install small baffle with no perforation at the rear of the upper tray support channel.

NOTICE

The perforated baffle is not to be used as a shelf to place media or other items. It is intended to eliminate splashing.

OPERATION CHECK

To check for proper operation of unit:

1. Close drain valve by turning handle clockwise.


WARNING

DO NOT OPEN DRAIN VALVE WHILE UNIT IS OPERATING. PREMATURE OPENING MAY RESULT IN SCALDING OF OPERATOR.


2. Fill chamber with 4 to 5 quarts (3.8 to 4.7 liters) of ordinary tap water. DO NOT USE DISTILLED OR DEIONIZED WATER.
3. Close chamber door.
4. Set exhaust selector to INSTRUMENTS AND PACKS (fast exhaust) or LIQUIDS (slow exhaust).
5. Set timer to 15 minutes. Cycle will go to completion automatically.

OPTIONAL COLD WATER CONDENSER

ITEM	PART No.	DESCRIPTION	QUANTITY
1	95-2119	Steam condensing unit	1
2	95-2219	Thermostat Box Assembly	1
3	95-0086	Exhaust line	1
5	15-7057	Copper tubing 3/8 OD, 22.25"	1
6	10-1775	Rd. Hd. Mach. Screw, 1/4-20	2
7	10-2500	Lockwasher, 1/4	2
8	10-2308	Hex Nut, 1/4-20	2
9	10-1812	Rd. Hd. Mach. Screw, 10-32	2
10	10-2505	Lockwasher, 10	2
11	10-2340	Hex Nut, 10-32	2
12	95-4009	Front Template (7" Lg)	1
13	95-4010	Back Template (11" Lg)	1




SERVICE CONNECTIONS REQUIRED	
A	1/2" IPS Cold Water Connection
B	1" IPS Drain Connection (See Note 3)
C	115V Elec. Connection 7/8 Ø knockout (cond. unit)
D	Electrical Connection




PAN SUPPORTS AND BAFFLES



1. Locate the mounting studs on the inside of the chamber. There are two rack mounting studs on each side.



2. Taking one pan support and positioning rack so that the pan stop is facing the rear of the unit and the wires are facing toward the center of the unit. The pan stop is a piece of sheet metal welded to the rack with a 65° bend.
3. Begin to hang the pan support by placing the rear key-way slot onto the rear mounting stud and slide the rack until the slot sits on the mounting stud. When this is done correctly the front mounting stud will be in position to place the front key-way slot. Slide the rack down into its correct position.



4. After installing one pan support rack correctly, you can install the upper baffle. Position the baffle so that the 45° bend is facing up towards the front of the unit (see "Figure 5"). Slide the mounting tab onto the flat bend on the pan stop bracket. The baffle should now stay in place by itself, but in a tilted state (see "Figure 6").
5. Position the second pan support rack into the cavity and slide the other mounting tab onto the rack flat bend while the pan support rack is not on the mounting studs. Hang the pan support by placing the rear key-way slot onto the rear mounting stud and slide the rack until the slot sits on the mounting stud. When this is done correctly the front mounting stud will be in position to place the front key-way slot. Slide the rack down into its correct position.




INSTALLATION

STERILMATIC - DIGITAL ELECTRIC STERILIZER

- Place the Perforated Water Baffle so that it sits on the bottom of the inside of the sterilizer chamber.

CAUTION

Do not manually refill hot, dry reservoir with cold water. Damage to the interior may result if cold water is added to a hot, dry reservoir.



Optional Stands

STERILMATIC OPEN STAND

Market Forge Sterilmatic Stand can be supplemented with an Optional Stand for utility use where maximum compactness is desired.

The sturdy, stainless steel unit is equipped with adjustable leg extensions which allow the unit to be installed and leveled over existing contours in the floor.

The open design lends itself to maximum sanitary conditions because of the ease with which periodic cleaning can be done.

Though simple in design and appearance, the sterilmatic stand is the ideal arrangement for mounting in that it allows secondary air to circulate.

STERILMATIC OPEN STAND WITH CONDENSER

Market Forge can provide the open stand with an optional steam condenser system for use where steam exhaustion into the room is undesirable.

The condenser is automatically controlled by the thermostat. The normal factory thermostat setting is 130°F (54°C). The open under-shelf of the stand gives added utility providing a handy tabouret for utensils and access for drainage of water from the sterilizing chamber.

PART No.	50 Hz	DESCRIPTION
10-4653	10-4653	Thermostat
10-4035	10-7074	3/8" Solenoid
10-5731	10-5731	1/2" Water Stop Valve
95-2106	95-2106	Water Injection Assy.
95-1680	95-1680	Shelf

GENERAL OPERATING INSTRUCTIONS

IMPORTANT

Make sure the drain valve is closed. Fill bottom of the sterilizer chamber with approximately six quarts of water or just below ledge at bottom of door opening. (If water supply is known to be hard or corrosive, a source of treated water should be used.) **DO NOT USE DISTILLED OR deionized WATER.**

1. **LOAD STERILIZER:** Use proper sterilizer loading procedures when placing materials in sterilizer chamber. All solid containers or instruments must be placed so that water or air will not be trapped in them.
2. **CLOSE DOOR:** Grasp handle, and holding it in vertical position, pull door down until bottom of door rests in the bottom of door opening. Then rotate handle forward, engaging the lower curved portion under the horizontal rod in the casting at the bottom of the door opening. Push handle all the way down and back until door is locked securely in position.
3. **DETERMINE CORRECT STERILIZATION TIMES:** Refer to the table on the following page for minimum sterilization times table.

NOTICE

In no case should the timer be set to less than 15 minutes. Sterilization will not be accomplished in less than 15 minutes exposure time.

4. When the sterilizer chamber reaches the selected temperature, the timed Heating/Sterilization cycle will begin. When the Heating/Sterilization cycle is completed, the electric supply to the heating elements will be opened (shut off) automatically. When the chamber pressure reaches 0 (zero) the door may be opened.

NOTICE

Before opening the chamber door be sure to have the Control Panel Flip Cover in the 'DOWN' position. This protects the LCD screen from coming into contact with too much steam. At this point you may release the handle and let go to avoid possible contact with the remaining escaping steam. When opening the door allow a few seconds for steam to escape from the chamber before opening completely.

NOTICE

For more detailed Operating Instructions please refer to "Detailed Operating Instructions" in this manual.

5. To assist in drying racks, release door handle after pressure has been attained at start of cycle. Pressure in chamber will keep door closed. The use of a wire basket will provide better drying for dressings. At end of sterilizing cycle, release door handle and open slightly. Do not lift door to open position. This will allow steam and moisture to escape. Allow door to remain in this position for 15 to 20 minutes before removing load. Small packs can be dried successfully with this procedure. We do not recommend the sterilization of large packs, such as linens. Be sure condensate baffles are in position in the chamber.
6. Remove load and check water level for next operation.

STERILIZATION GUIDE

- **PACKS** (Linens, gloves, etc.): Use wire basket to facilitate drying. Be sure condensate baffles are in place. Place packs on edge and arrange load in chamber, so that only minimal resistant to passage of steam through the load will exist.

NOTICE

Place gloves in upper two-thirds of chamber.

- **JARS, CANISTERS** (etc.): Place containers on side to allow for displacement of air and complete contact of steam to surfaces. Drying is also facilitated.
- **PETRI DISHES, PIPETTES, DESICCATORS** (etc.): Should be inverted.
- **UTENSILS, TREATMENT TRAYS**: Placed on edges to facilitate drying.
- **INSTRUMENT SETS**: Place instruments set in trays having mesh or perforated bottoms. Place trays flat on shelves.
- **COMBINING FABRICS & HARD GOODS**: Place hard goods on lowest shelves.
- **PLASTIC UTENSILS**: DO NOT stack or nest plastic items.
- **LIQUIDS**: Sterilize liquids separately from other supplies or materials. Set vent to slow.
- **SMALL ITEMS**: Sterilize small items in baskets, or trays.

IMPORTANT

IF THE EQUIPMENT IS USED IN A MANNER NOT SPECIFIED BY THE MANUFACTURER, THE PROTECTION PROVIDED BY THE EQUIPMENT MAY BE IMPAIRED.

MINIMUM STERILIZATION TIMES

TIME (MINUTES)	ARTICLES
15	<ul style="list-style-type: none"> • Glassware, empty, inverted. • Instruments, metal in covered or open tray, padded or unpadded. • Needles, unwrapped. • Pipettes, blood diluting, serological, volumetric, etc • Tubing glass (6mm), (10mm) inverted
20	<ul style="list-style-type: none"> • Flasked solutions 75-250 ml. • Instruments, metal combined with other materials in covered and/or padded tray. • Instruments wrapped in double thickness muslin. • Rubber gloves, catheters, drains, tubing, etc. Unwrapped or wrapped in muslin or paper.
30	<ul style="list-style-type: none"> • Brushes in dispensers, in cans of individually wrapped. • Dressings, wrapped in paper or muslin, small packs only. • Flasked solutions 500-1000 ml. • Syringes, unassembled, individually packaged in muslin or paper. • Needles, luer, individually packaged in glass tubes or paper.
45	<ul style="list-style-type: none"> • Flasked solutions 1500-2000 ml.

DIGITAL CONTROL PANEL

Figure 8




Figure 9

Before operating this unit, be sure to have read the Owner’s manual for proper setup, service connections and installation. In addition, the Owner’s manual will cover sterilization recommendations, daily cleaning procedures and parts lists.

A fingertip or stylus can be safely used to navigate the touch screens. A detailed explanation of the functions of each with step-by-step procedure of operations is given below.

NOTICE

1. Never squirt water on the digital panel.
2. Do not operate the console keys using a sharp- pointed object such as a propelling pencil or needle. Doing so might damage the console.

ITEM	DESCRIPTION
1	Control Panel Power - ON (1), OFF (0)
2	USB Port
3	Printer Input
4	Reset Button
5	Low Water Indicator Light
6	Digital LCD Display
7	Power Key
8	Settings Key
9	LCD Cleaning Screen Lock Key
10	Save Key
11	Forward Key
12	Cancel Key
13	Print Key
14	View Key
15	USB Key

The Digital Control panel is made up of Operational Switches, Digital LCD Touch Screen and external Data retrieve ports for the Printer and USB.

For the sterilization process, the operator can either choose ‘Manual’ or ‘Presets’ modes. In both cases, the Digital LCD Display will show the actual temperature within the sterilizer chamber, the time remaining for the sterilizer cycle and the chamber pressure in digital form. It also displays the current process state, whether Preheating, Sterilizing or Venting.

The three operating parameters that the operator can set are sterilizing Temperature, Time, and Venting mode. The temperature is adjustable in the range between 225°F (107°C) to 250°F (121°C). The sterilizing time can be adjusted in minutes and seconds. The Fast-venting mode takes approximately 3 minutes, whereas the Slow- venting mode takes approximately 11 minutes.

After initiating the Manual Sterilizing Cycle, the operator will be able to set three inputs. Temp, Time, and Vent. Or the operator can save commonly used values for these parameters into ‘Presets’ keys. The parameter values stored in the Preset keys can easily be selected before running a sterilization cycle. This avoids the operator having to set all three parameters every time the operator runs the unit.

SETUP

- Make sure the unit is connected to its electrical source.
- Make sure the unit has recommended amount of water in its chamber.
- Make sure the unit drain is closed. If a ventilation hood is required. Make sure unit is placed accordingly under the hood.
- Be sure that the unit's door is closed and locked before operating.

Once the unit has power, push the black Control Panel Power switch to ON (1). This will bring power to the control panel and illuminate the Digital LCD Display screen. The unit is now ready to use.

STEP BY STEP PROCEDURE

MAIN SCREEN

By selecting the Settings Key (8) the operator has access to Setup and Service selection options. By selecting the Setup, the operator can change the parameters like units (°F/PSI to °C/KPa), Time & Date and Buzzer volume (Low, Medium & High). The Service option is for manufacturer's use only. In case of any technical issues, contact the manufacturer service department.

By pressing and holding the LCD Cleaning Screen Lock Key (9) for 2 seconds, the operator can lock the digital screen for the next 20 seconds. Which permits the operator to clean the LCD screen without touching or modifying any functions during operation.

By selecting the Power Key (7) at the center of the display, the operator can access the following options: Presets, Manual.

The operator can also provide USB/printer requests on the same screen and go back to the main menu by touching the Power Key. Detailed explanations of both the options and set up are given in this section.

Figure 9





Figure 10



Figure 11




Manual Menu

The Manual setup screen provides options to set desired Times, Temperatures and Slow or Fast vent. The values displaying on the opening screen are the default values.

SETTING TIME, TEMP & VENT

1. At the starting of the cycle, the displayed Time and Temp are at default values. The operator can adjust these to desired sterilizing Time and Temp by pressing and holding ' + ' & ' - '.
2. In the 'Idle State,' the Vent mode will be displayed as either Fast or Slow, whichever is active. To change the Vent mode, simply press the Slow or Fast option.
3. After adjusting to desired values, using the Forward Key (11) operator can directly start the process.



PRESETS MENU

In the Presets menu, the values of Time, Temp and Vent speed stored are standard pre-defined values. Compare to the Manual setting Presets option will save a lot of time during sterilization operation.


SETTING TIME, TEMP & VENT

Programming of Time, Temp and Vent in Presets mode follows the same procedure of Manual programming.

However, for setting and saving a new program, the operator must go to the Manual menu.

In the Manual Menu Screen,

1. At the starting of the cycle, the displayed Time and Temp are at default values. The operator can adjust these to desired sterilizing Time and Temp by pressing and holding ' + ' & ' - '.
2. In the 'Idle State,' the Vent mode will be displayed as either Fast or Slow, whichever is active. To change the Vent mode, simply press the Slow or Fast option.
3. After adjusting to desired values, using Save key (10), save the values under new name. Return to the Presets Screen to view saved data.



STARTING THE UNIT

1. Once you have set your desired Temp, Time and Vent type you can now proceed in starting the unit.
2. For starting in the manual mode of operation, select the Forward Key (11) that appears on the Sterilization Parameters screen, or for starting using presets in the Presets Selection screen, by selecting the pre-saved item.

STOPPING THE UNIT

1. Once the unit is started it will run through all three cycles automatically (cycles are described in this section) and finish in the Done state (also described in this section).
2. If at any time throughout the cycle you must stop or cancel the cycle just press the Cancel Key (12).
3. After stopping the unit in this manner, the unit will still have to go through the Vent cycle to release the pressure within the chamber before you can open the door.

STERILIZING CYCLES

1. **Pre-Heating Cycle:** Once the manual or presents key is selected, the sterilizer is heating up to the sterilization temperature set point. While heating, the set point time is displayed and will not change, but the chamber's temperature and pressure gradually increases.
2. **Sterilizing Cycle:** When the unit reaches its set-point temperature the unit will enter the sterilization cycle. At this point the unit will display and count down the sterilization time. At this point the unit will display and count down the sterilization time.
3. **Venting Cycle:** When the Sterilization Cycle time reaches zero, the set Venting Cycle begins. The Fast Vent is programmed to vent for 3 minutes through the solenoid valve. The Slow Vent is programmed to vent for 11 minutes through a bleeder orifice.
4. **DONE state** – After the Venting Cycle, the unit will display Done, and a beeper will sound. The beeper sequence will be ON for 1 second and OFF for 9 seconds. This sequence will repeat until any key on the control panel is pressed or 3 minutes have passed.
 - Once the system displays a Done state, it is safe to open the chamber door. Before opening the Chamber door, be sure to have the Control Panel Flip Cover in the 'Down' position. This protects the LCD screen from getting into contact with too much- escaping steam from the chamber. When the door is left in the covered position, the operator may notice a rise in the temperature and pressure on the LCD screen. The operator should open the chamber door to allow the chamber to cool.
 - The operator can view complete cycle information in Done screen using View Key (14). Also, the operator can retrieve cycle data by either pressing USB Key (15) or Printer Key (13) on the same screen.

NOTICE

There is never any harm in releasing the chamber door latch. If there is pressure inside the chamber, the door cannot open due to its design which does not allow the door to open under pressure.

Figure 14




Figure 15





Figure 16



Figure 17



PRINTING/USB DATA LOGGER

By connecting the USB Flash Drive to the USB port located on the panel (Refer Figure: 8) (2), the operator will be able to retrieve the Time, Temperature and Pressure during all phases of the entire sterilization cycle in an excel file format. All three of these parameters are recorded at one-minute intervals. The data logging clearly labels the output with headers separating the three cycles, Heating, Sterilizing and Venting.

Each run is assigned a unique ID number that is displayed at the top of the data sheet. Each unique ID number is made up of two sections separated by a dash (-) (Example; #LAB-00 005). The section (LAB) never changes. It is a unique number that identifies the sterilizer unit itself. The last five digits are incremented by one for every time the unit is run.

Also, if the operator wishes to include any other information like the Operator Name or Comments on the retrieved data, the operator can easily edit the excel sheet.

NOTICE

The sterilizer has a connection port for an optional external printer (not available from the manufacturer) at the control panel (Refer Figure: 9) (3). The external printer needs to have an RS232 connection and be compatible with ESC/POS command.

RETRIEVING PREVIOUSLY RUN CYCLES.

The data logger will store approximately 20 complete sterilization cycles. This number may vary slightly depending on how long each cycle runs. For example, the data logger can store more cycles set at 30 minutes versus cycles set at 60 minutes. As stated before, each complete cycle is assigned a unique ID number. Thus, the operator can view full sterilization cycle data on the LCD screen by selecting the desired Id. In addition, on the same screen, the operator will be able to give print commands or save data to a USB drive (Refer: Figure 18)

Figure 18

File Exported on 07/21/21 @ 11:55:50

LAB-00005

Date: 07/19/21

Time: 10:45AM

Preset: 250FAST

Set Time: 10:00

Set Temp: 250°F

Heating

0:00	198.0°F	0.0PSI
1:00	230.7°F	6.4PSI
2:00	247.3°F	13.8PSI
3:00	249.7°F	15.0PSI
3:06	250.0°F	15.1PSI

Sterilizing

0:00	250.0°F	15.1PSI
1:00	250.8°F	15.6PSI
2:00	250.3°F	15.3PSI
3:00	249.8°F	15.0PSI
4:00	249.8°F	15.0PSI
5:00	250.0°F	15.1PSI
6:00	250.1°F	15.2PSI
7:00	250.0°F	15.1PSI
8:00	250.0°F	15.1PSI
9:00	250.0°F	15.1PSI
10:00	250.0°F	15.1PSI

Venting

0:00	250.0°F	15.1PSI
1:00	214.6°F	0.8PSI
2:00	216.4°F	1.3PSI
3:00	217.0°F	1.5PSI

LOW WATER RESET

If the water inside the chamber is allowed to run dry it will trigger the Low Water Cut-off. At this point the unit will;

- Shut down all three heating elements
- Light the red Low Water indicator light.
- Start the Venting cycle
- The Data Logger will record the error code

NOTICE

The sterilizer's LCD display may show the temperature and pressure rising slightly during this venting period but that is normal.)

- Once the unit completes the VENTING cycle the screen will flash 'DONE' and display 'SERVICE' and the Beeper will sound.

STEPS TO RESET THE UNIT

1. The LCD display now shows 'ERR 06' and beeper changes to a constant tone.

NOTICE

An explanation of Error Codes is shown below in the Appendix.

RECOMMENDED STEPS BUT NOT REQUIRED:

1. Shut off 'ON/OFF' power switch to controller.
2. Wait until unit cools down. Opening the unit door will help the unit to cool down quicker.
3. When power is restored the buzzer tone will continue and 'ERR 06' will be displayed again.
4. Open the door and add water to the chamber.
5. Once the chamber has cooled enough you can press the Low Water Cutoff RESET button. Listen for the 'click' sound.
If constant tone continues press any menu key to turn off.
6. The unit is now back in the 'Idle State' and ready to be run again.

MAINTENANCE

CLEANING

DAILY CLEANING PROCEDURE (AT THE END OF EACH DAY)

When sterilizing culture mediums that generate sulfide gas or chlorine gas, the inside of the chamber must be cleaned and rinsed thoroughly without fail.

1. Remove bottom splash baffle.

IMPORTANT

STERILIZING CHAMBER MUST BE CLEANED AND DRAINED DAILY USING THE FOLLOWING PROCEDURE. WASH WETTED PORTION OF THE CYLINDER THOROUGHLY BY ADDING A MILD DETERGENT TO WATER IN CYLINDER.

2. If a soft cloth or brush is used with the detergent and does not completely remove the surface film, a nylon soap pad should be used. After washing thoroughly rinse with clean water. Dry cylinder* and leave door open overnight.

* The Sterilmatic cylinder is constructed of corrosion resistant Alclad aluminum alloy. The protective properties of this material afforded to the interior portion of the cylinder which is exposed to water may be destroyed by allowing a film to form. Such a film can be caused by salts or other contaminants in the water. Corrosion may also occur if water is not drained daily.

WEEKLY CLEANING

In addition to the daily cleaning it is necessary to clean the air intakes on a weekly basis. Air intakes provide necessary cooling air to the internal components. They are generally located on the rear and sides of the equipment.

Warranty

STERILIZER (AUTOCLAVE) WARRANTY

MODELS: STM-ED, STM-EDX*

We warrant to the original purchaser that the sterilizers manufactured by Market Forge will be free from defects in material and factory workmanship if properly installed and operated under normal conditions. Within one year from date of original installation, or within 15 months from date of shipment from factory, whichever is sooner, we will repair or replace that part of any such machine that becomes defective at no cost to the customer.

This warranty is effective for One (1) Year Parts and 90 Days Labor, Travel and Mileage.

This warranty does not apply to damage resulting from use of hard or corrosive water, from failure to drain and dry cylinder daily or from inadequate cleaning procedures. Nor does it cover any part or assembly, which has been subjected to accident, alteration, or is from a machine where the serial number has been removed or altered. Normal service adjustments are not covered by this warranty.

Any defect during the warranty period shall be brought to the attention of a factory authorized service agency or the dealer from whom the equipment was purchased. He will be authorized to furnish or arrange for repairs or replacements within the terms of the warranty.

NOTICE

This warranty only applies to the USA and Canada. Elsewhere, warranty covers parts only for one year as described above.

* Export Model.



MARKET FORGE

Telephone: (919) 762-1000

www.crownsteamgroup.com

NOTES

STM-ED & STM-EDX Sterilmatic Digital Electric Sterilizer



A product with the Market Forge name incorporates the best in durability and low maintenance. We all recognize, however, that replacement parts and occasional professional service may be necessary to extend the useful life of this appliance. When service is needed, contact a Market Forge Authorized Service Agency, or your dealer. To avoid confusion, always refer to the model number, serial number, and type of your appliance.



CROWN FOOD SERVICE EQUIPMENT

 **MIDDLEBY** A Middleby Company

70 Oakdale Road, Downsview (Toronto) Ontario, Canada, M3N 1V9
Telephone: 919-762-1000 www.crownsteamgroup.com