

STERISURE™

Precision Controlled Dry Heat Sterilizers

Models 2100 and 3100



Model 2100

Made in USA



Optional printer

SteriSURE Dry Heat Sterilizers are microprocessor controlled units for all metal and steam sensitive instruments.

- ✓ Audible and visual indication of warm-up and completion of sterilization cycle
- ✓ Paperless data storage
- ✓ Error message will display if cycle is interrupted
- ✓ Temperature can be displayed in Celsius or Fahrenheit
- ✓ Automatic shutoff (adjustable from 2-4 hours after cycle is complete)
- ✓ Maintains temperature between cycles
- ✓ Solid State control of heaters

One Touch Automatic Controls and Digital Display

SteriSURE Dry Heat Sterilizers automatically start the sterilizer cycle and set the temperature at 320° F (160° C) and time at 60 minutes. The time and temperature are displayed throughout the cycle allowing the operator updated sterilization information.



Advantages of Dry Heat Sterilization

- SteriSURE Dry Heat Sterilizers kill bacteria and spores without moisture.
- Does not dull, pit, or rust sharp instruments.
- Dry heat sterilization does not require a dry cycle that adds to the sterilization cycle time.
- Economical, little to no maintenance, low operating costs and requires no distilled water.

cpac
Equipment Inc.

STERISURE™

2364 Leicester Road • PO Box 175 • Leicester, NY 14481
Phone: 1-800-828-6011 • (585) 382-3223 • Fax (585) 382-9481

www.sterisure.com

Precision Controlled Dry Heat Sterilizers



Made in USA

Optional printer



Model 3100



Printer Output/USB Port

Both models of the SteriSURE Dry Heat Sterilizer have a USB port to provide a permanent record of the sterilization cycle.

Maintenance Free

SteriSURE Dry Heat Sterilizers require no routine cleaning. The high quality stainless steel construction maintains a great looking unit for years.

Specifications

- Sterilization Time: 60 minutes
- Unit Size:
 - Model 2100 15.25"W x 11.75"H x 11"D, 387 mm x 298 mm x 279 mm
 - Model 3100 17.75"W x 12.75"H x 13"D, 450 mm x 324 mm x 330 mm
- Chamber Size:
 - Model 2100 11.25" x 4.25" x 7", 286 x 108 x 178 mm
 - Model 3100 13.75" x 6.25" x 8", 349 x 159 x 203 mm
- Weight:
 - Model 2100 17 lbs (8 Kg)
 - Model 3100 27 lbs (13 Kg)
- Dimensional Shipping Weight/Dimensions:
 - Model 2100 20 lbs, 20" x 16" x 14"
 - Model 3100 39 lbs, 22" x 18" x 16"
- Tray Size:
 - Model 2100 (2 trays) 10"W x 6.25"D, 254 mm x 159 mm
 - Model 3100 (3 trays) 11.75"W x 7.25"D, 298 mm x 184 mm
- Electrical: 110V-120V, 60Hz, Model 2100, 500 watts
220V-240V, 50Hz, Model 3100, 750 watts
- Sterilization Temperature: 320° Fahrenheit (160° Celsius)
- Warranty: One year (parts and labor)
- Models: 2100, 3100
- Approvals: ETL (ETL logo) FDA 510(k) K094026
IEC Safety Standard 61010-2-040
Validated to AAMI Standards ST-40 and ST-50
UL 61010-1, CAN/CSA C22.2 No. 61010-1



cpac
Equipment Inc.

2364 Leicester Road
PO Box 175
Leicester, NY 14481
Phone: 1-800-828-6011
+1 (585) 382-3223
Fax +1 (585) 382-9481
www.sterisure.com

Authorized Distributor:

