

Sponsor: Pierre Dizier GVS Filter Technology UK Vickers Industrial Estate Mellishaw Lane Morecambe, Lancashire, LA3 3EN UNITED KINGDOM

Bacterial Filtration Efficiency (BFE) and Differential Pressure (Delta P) Final Report

F31000/N31000 Test Article:

Lot #OL 20/00551

Purchase Order: FO-2020/01722/21 Study Number: 1299170-S01

Study Received Date: 13 May 2020

Testing Facility: Nelson Laboratories, LLC 6280 S. Redwood Rd.

Salt Lake City, UT 84123 U.S.A.

Standard Test Protocol (STP) Number: STP0004 Rev 18 Test Procedure(s):

Deviation(s): None

Summary: The BFE test is performed to determine the filtration efficiency of test articles by comparing the bacterial counts upstream of the test article to the bacterial counts downstream. A suspension of Staphylococcus aureus was aerosolized using a nebulizer and delivered to the test article at a constant flow rate and fixed air pressure. The challenge delivery was maintained at 1.7 - 3.0 x 10³ colony forming units (CFU) with a mean particle size (MPS) of $3.0 \pm 0.3 \mu m$. The aerosols were drawn through a sixstage, viable particle, Andersen sampler for collection. This test method complies with ASTM F2101-19 and EN 14683:2019, Annex B.

The Delta P test is performed to determine the breathability of test articles by measuring the differential air pressure on either side of the test article using a manometer, at a constant flow rate. The Delta P test complies with EN 14683:2019, Annex C and ASTM F2100-19.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

> Test Side: Inside

~9.1 cm² BFE Test Area:

BFE Flow Rate: 28.3 Liters per minute (L/min) Delta P Flow Rate: 8 Liters per minute (L/min)

Conditioning Parameters: $85 \pm 5\%$ relative humidity (RH) and 21 ± 5 °C for a minimum of 4 hours

Test Article Dimensions: ~230 mm x ~155 mm

1.7 x 10³ CFU Positive Control Average:

Negative Monitor Count: <1 CFU

> MPS: 3.2 µm





David Brown electronically approved for

James Luskin

05 Jun 2020 05:03 (+00:00) Study Completion Date and Time

Study Director

801-290-7500

nelsonlabs.com

sales@nelsonlabs.com

FRT0004-0001 Rev 22 Page 1 of 2



Results:

Test Article Number	Percent BFE (%)
1	>99.9 ^a
2	>99.9 ^a
3	>99.9
4	>99.9 ^a
5	>99.9 ^a

^a There were no detected colonies on any of the Andersen sampler plates for this test article.

Test Article Number	Delta P (mm H ₂ O/cm ²)	Delta P (Pa/cm ²)
1	12.6	123.3
2	11.7	114.9
3	10.3	101.2
4	12.9	126.1
5	11.8	115.8

The filtration efficiency percentages were calculated using the following equation:

$$\% BFE = \frac{C - T}{C} \times 100$$

C = Positive control average

T = Plate count total recovered downstream of the test article Note: The plate count total is available upon request



Sponsor:
Rob Davis
GVS Filter Technology UK
nfc house, vickers ind estate, mellishaw
morecambe, lancashire, LA3 3EN
UNITED KINGDOM

Synthetic Blood Penetration Resistance Final Report

Test Article: F30000/N31000 Facemasks

Purchase Order: FO-2020/02277/21 Study Number: 1311153-S01 Study Received Date: 17 Jun 2020

Testing Facility: Nelson Laboratories, LLC

6280 S. Redwood Rd.

Salt Lake City, UT 84123 U.S.A.

Test Procedure(s): Standard Test Protocol (STP) Number: STP0012 Rev 09

Deviation(s): None

Summary: This procedure was performed to evaluate surgical facemasks and other types of protective clothing materials designed to protect against fluid penetration. The purpose of this procedure is to simulate an arterial spray and evaluate the effectiveness of the test article in protecting the user from possible exposure to blood and other body fluids. The distance from the target area surface to the tip of the cannula is 30.5 cm. A test volume of 2 mL of synthetic blood was employed using the targeting plate method.

This test method was designed to comply with ASTM F1862 and ISO 22609 (as referenced in EN 14683:2019 and AS4381:2015) with the following exception: ISO 22609 requires testing to be performed in an environment with a temperature of $21 \pm 5^{\circ}$ C and a relative humidity of $85 \pm 10\%$. Instead, testing was performed at ambient conditions within one minute of removal from the environmental chamber held at those parameters.

All test method acceptance criteria were met. Testing was performed in compliance with US FDA good manufacturing practice (GMP) regulations 21 CFR Parts 210, 211 and 820.

Number of Test Articles Tested: 32 Number of Test Articles Passed: 32

Test Side: Outside

Pre-Conditioning: Minimum of 4 hours at 21 ± 5°C and 85 ± 5% relative humidity (RH)

Test Conditions: 23.5°C and 21% RH

Results: Per ASTM F1862 and ISO 22609, an acceptable quality limit of 4.0% is met for a normal single sampling plan when ≥29 of 32 test articles show passing results.

Test Pressure: 160 mmHg (21.3 kPa)

<u> </u>	•
Test Article Number	Synthetic Blood Penetration
1-32	None Seen



Adam Brigham electronically approved for

22 Jul 2020 15:24 (+00:00)

Study Director

James Luskin

Study Completion Date and Time

801-290-7500 | nelsonlabs.com | sales@nelsonlabs.com

FRT0012-0002 Rev 13 Page 1 of 1