

STERILIZATION CENTERS

MAXIMUM FUNCTIONALITY FOR INSTRUMENT REPROCESSING

SEFFICIENT.



品 SYSTEMATIZED.

ॐ ATTRACTIVE.

RESILIENT.

THE STERILINE™ STERILE CENTER IS DESIGNED TO GIVE MAXIMUM FUNCTIONALITY DURING THE CRITICAL STAGE OF INSTRUMENT REPROCESSING - FROM CLEAN-TO-STERILE.



DEPOSIT WASTE

TEMPORARY WASTE STORAGE

Removable waste container. Push-pull door allows hands-free operation.

CLEAN PROCESS

ULTRASONIC CLEANING / MANUAL CLEANING AREA The DURAsonic cleaners can

The DURAsonic cleaners can be elevated off the work top for additional work space.

WASHER DISINFECTOR

HIGH CAPACITY CLEANING Space allows for all brands, including Miele, Hydrim, or Melag brand washers.

BAGGING + PACKAGING

MOUNTED DISPENSER
Keeps tubing off counter.
Storage cabinets keep autoclave wrap or pouches/tubing organized.

STERILIZATION

AUTOCLAVE
Sterilizer can be either inserted into the cabinet below – or placed on work top.

STORAGE

STERILE STORAGE

Base and wall units designed to neatly store sterile instruments.







Suspend medical waste container.



Suspend autoclave tubing dispenser.







Push-pull doors for hands free operation. Removable bins.




Increase worktop area when needed.




WALL UNITS



Hygiene center includes towel, mask and soap dispensers.





two hinged doors (two shelves inside)

Wall unit with

WORK TOP ORGANIZER



WORK TOPS





- - 3 | Stratified laminate worktop

1 | Corian® worktop

2 | Stainless steel top

4 | Stratified "Fenix" laminate worktop







Work top pop-up power station and data recepticles.

BASE UNITS





Choose between closed base or modern aluminum foot base.


15

STERILINE COMPOSITIONS

CUSTOMIZABLE TO SUIT YOUR SPACE.

STERILINE 2715









CALL TOLL FREE 877-561-0500 VISIT DuralineSystems.com EMAIL sales@duralinesystems.com

STERILINE 3610





STERILINE 2110



STERILINE 2210

