

A TRUE ADVANCEMENT IN STERILIZATION HAS ARRIVED!

COX RAPIDHEAT

CPAC
ENVIRONMENTAL SOLUTIONS

HIGH-VELOCITY HOT AIR (HVHA) STERILIZER

**Rapidly Being Proven as the Most Cost-Effective,
Environmentally-Friendly Sterilization Method in the World!**

EXAMINE THE HVHA ADVANTAGE:

- ▶ Complete FDA-Cleared 12-Minute Wrapped Cycle
- No Shortcuts for Immediate-Use
- ▶ Uses 84% Less Energy Than Steam
- Verified Through Independent Research
- ▶ Does Not Corrode, Pit, or Dull Instruments
- Extends Life of your Instruments
- ▶ Does Not Use Water, Steam or Chemicals
- An Environmentally Sustainable Process
- ▶ Little or No Maintenance & Repair
- Two Year Savings Can Pay for HVHA



MADE IN THE
USA



CoxRAPIDHEAT
HIGH-VELOCITY HOT AIR STERILIZERS

**NEW THREE-YEAR
PARTS & LABOR WARRANTY**

HVHA Sterilizers save money, time and resources – It's the most efficient sterilizer solution for any practice!

Contact your favorite dealer, call 800-828-6011 or visit www.CPAC.com.



CoxRapidHeat™

HIGH-VELOCITY HOT AIR STERILIZERS

Over 25 Years of Professional Use

Specifications:

Unit Size (OD)	17" W x 17"D x 15"H, 432 mm x 432 mm x 381 mm
Chamber Size (ID)	9.625" W x 7.75" D x 4" H, 244 mm x 197 mm x 102 mm
Weight Unit	50 lbs (23 Kg)
Dimensional Shipping	61 lbs, 24" x 22" x 19"
Instrument Tray Size (OD)	9.5" W x 7.75" D x 1.875" H, 241 mm x 197 mm x 48 mm
Electrical	110-120V, 60Hz, 1100 Watts warm up, 300 Watts operating 220-240V, 50Hz, 1100 Watts warm up, 300 Watts operating
Sterilization Temperature	375° Fahrenheit (190° Celsius)
Warranty	Three-Year (parts and labor)
Model	6000-16
Approvals	ETL FDA 510(k) K881371 Health Can Lic #93771 UL 61010-1 CAN/CSA C22.2 No.61010-1, CSA C22.2 No. 61010-2-040-07



Sterilization Times

Cycle 1: Six (6) minutes for unwrapped instruments


Cycle 2: Eight (8) minutes for hand pieces and medical drills

Cycle 3: 12 minutes for packaged instruments

Included with Sterilizer



Small Cassette Rack



Instrument Pouch Rack



Tray Removal Tool



Cooling Rack



Instrument/Rack Tray



2364 Leicester Road • PO Box 175
Leicester, NY 14481
(800) 828-6011 • (585) 382-3223
www.CPAC.com

*For a full complement of sterilization supplies
(pouches, spore test kits and indicator strips)
visit www.CPAC.com*

